



NON-SLIP DESIGN

The next generation in
airway balloon dilation



aeris is the most advanced balloon catheter for airway dilation offering superior design to address safety concerns. Conceived by an airway surgeon, no other balloon catheter specifically addresses the issues of slippage in airway stenosis.

Truly a safe, efficient answer for treatment of airway stenosis.



atraumatic access



safe, secure placement



controlled radial dilation



consistent dilation of strictures

Non-slip Design

Provides safe, secure placement before dilating, reducing the risk of slippage.

Single Lumen Catheter

Leads to rapid inflation and deflation and maximum efficiency in a time sensitive procedure.

Wide Range of Sizes

Allows for optimal fit to accommodate both adult and pediatric airway diameters.

Cost Effective

Meets the needs of today's healthcare economics.

Catalog Number	Units Per Box	Inflated Balloon Diameter (mm)	Balloon Length (mm)	Recommended Inflation Pressure (atm)
KG0530	5	5	30	17
KG0730	5	7	30	17
KG0830	5	8	30	17
KG0930	5	9	30	17
KG1030	5	10	30	17
KG1240	5	12	40	10
KG1440	5	14	40	10
KG1640	5	16	40	10
QL2530	Disposable Inflation Device	One each per box		



EMAIL: sales@bryanmedical.net

WEB: www.bryanmedical.net

PHONE: 513-272-1600

Patent. US7,771,446B2