



Designed for the treatment of swallowing impairment caused by upper esophageal sphincter dysfunction.

- Simple transnasal delivery system
- Novel design tailored to the anatomic configuration of the UES
- Contoured to the specific physiology of the UES, allowing dilation to nearly double the inflation diameter of conventional devices





Infinity[®] has a 38 mm maximum outer diameter and **nearly double** the inflation dimension of conventional devices.

INFLATED INFLATED PRODUCT **BALLOON** INFLATION BALLOON NUMBER PRESSURE DIAMETER LENGTH atm mm mm **INFINITY 1000 INF1023** 5 atm 23 mm 55 mm THE FARWELL **INFINITY 3000** 32 mm **INF3032** 4 atm 55 mm THE MERATI **INFINITY 5000** 4 atm 38 mm 55 mm **INF5038** THE POSTMA

SPECIFICATIONS



Available exclusively from Bryan Medical 6100 Wooster Pike Cincinnati, OH 45227 513.272.1600 bryanmedical.net

HOPE MEDICAL

hopemedicalinc.com