

Order Form for TRACOE *comfort* / *comfort XL*
Product variants*

This is a digital fillable PDF

Your Information

Customer	
Name and contact details	
Patient ID / No. (anonymised)	

Order

REF + size	Number of units
Has been ordered before / Repeat order	
Yes ☐	Order number.:
No ☐	

Outer cannula (see explanations on pages 2 and 3)

Length variant may deviate from the standard length by max. +/- 10 mm).

By altering the length, the ID / OD of the tube may be smaller or larger than stated in the catalogue, due to its conical design. Please note that only the TL is given on the outer packaging.

AB: TL:

Inner cannula(s)


☐ Standard (as in catalogue) ☐ Other

Number of additional inner cannulas:

Inner cannula	Fenestrated	Un-fenestrated	15 mm connector	22 mm Adapter	Base	Speaking valve		Additional oxygen supply port
						Type A	Type B	
								
IC-1	☐	☐	☐	☐	☐	☐	☐	☐
IC-2	☐	☐	☐	☐	☐	☐	☐	☐
IC-3	☐	☐	☐	☐	☐	☐	☐	☐
IC-4	☐	☐	☐	☐	☐	☐	☐	☐

Total number of inner cannulas:

Bending angle and hook position



Standard bending angle	Alternative bending angle	Hook position	
☐ 90°	☐ 100°	☐ bottom	☐ top
	☐ 110°	☐ standard	
	☐ 120°		

Fenestration

Please select the desired fenestration. Please state if you require the standard distance to the start of the fenestration (L1) or whether you require a different distance. The following applies to all variants: the minimum distance for L1 is 10 mm for all can-nula sizes. You can increase the standard L1 by a maximum of 15 mm to match your desired tube size. This increase can be done in 2.5 mm steps.

☐ **Standard fenestration:**
☐ **Alternative fenestration:**
☐ **No fenestration**

☐ **L1 standard**
☐ **L1 variable:** _____ mm
 ☐ **L1 standard**
☐ **L1 variable:** _____ mm

Standard fenestration of TRACOE *comfort*:

Size	03	04	05	06	07	08	09	10	11	12	13	14
Hole pattern												
Hole dia-meter Ø	2.0	2.0	2.5	2.5	3.0	3.0	3.5	4.0	4.0	5.0	5.0	5.0
Length L2	7	7	14	14	15	15	17	18	18	21	21	21
Width L3	5	5	7	7	8	8	10	11	11	14	14	14
L1	10	10	12,5	15	17,5	20	22,5	25	27,5	27,5	27,5	27,5
TL	42	50	52	54	57	61	65	74	76	78	87	91


Standard fenestration of TRACOE *comfort XL*:

Size	03	04	05	06	07	08	09	10	11	12	13	14
Hole pattern	n.a.	n.a.										
Hole dia-meter Ø			2.5	2.5	3.0	3.0	3.5	4.0	4.0	5.0	5.0	5.0
Length L2			15	15	15	13	15	18	18	18	18	18
Width L3			7	7	8	8	10	11	11	14	14	14
L1			12,5	15	17,5	20	22,5	25	27,5	27,5	27,5	27,5
TL			85	87	93	96	100	104	108	112	116	122

Alternative fenestration of TRACOE *comfort* and TRACOE *comfort XL*:

Size	03	04	05	06	07	08	09	10	11	12	13	14
Hole pattern												
Hole dia-meter Ø	2.0	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5
Length L2	6,5	8	11,5	11,5	14,5	17,5	17,5	20,5	20,5	20,5	20,5	20,5
Width L3	4,5	5,5	6	6	6	6	9	9	9	9	9	9
L1	10	10	12,5	15	17,5	20	22,5	25	27,5	27,5	27,5	27,5

Explanation of the dimensions



- TL:** Length along the centre line from the start of the neck flange to the lower end of the cannula
- AB:** Length along the outer bend (greatest length)
- L1:** Distance between the neck flange and the first hole of the fenestration
- L2:** Length of the fenestration
- L3:** Width of the fenestration
- H:** Position of hook (top or bottom)
- Θ°** Bending angle

Comments

If you have any questions, please contact us by phone at (513) 272-1600

Please fill out the order form and send it to sales@bryanmedical.net or by fax to (513) 272-1610

Date _____ Signature and company stamp _____

