



TRACOE silcosoft[®] *Custom Program*

Neonates, Infants & Children

STERILE • IN STOCK • READY TO SHIP

When a standard size tracheostomy tube doesn't offer the right fit, turn to the TRACOE silcosoft[®] Custom Program.

- TRACOE silcosoft[®] has 65 custom tube lengths and 4 cuff styles designed to fit neonates, infants and children.
- Address patient needs sooner and reduce admission days with timely delivery of custom length TRACOE silcosoft[®] tracheostomy tubes.



TRACOE[®]

Standard Lengths (mm)

TO ORDER, SPECIFY: REF # + LENGTH (e.g. 360 - 2.5)

Standard, without H2O Cuff



Neonatal: REF 360

Pediatric: REF 370

| LENGTHS | ORDER # |
|----------|-----------|
| 2.5 x 30 | 360 - 2.5 |
| 2.5 x 38 | 370 - 2.5 |
| 3.0 x 32 | 360 - 3.0 |
| 3.0 x 39 | 370 - 3.0 |
| 3.5 x 34 | 360 - 3.5 |
| 3.5 x 40 | 370 - 3.5 |
| 4.0 x 36 | 360 - 4.0 |
| 4.0 x 41 | 370 - 4.0 |
| 4.5 x 42 | 370 - 4.5 |
| 5.0 x 44 | 370 - 5.0 |
| 5.5 x 46 | 370 - 5.5 |

Standard, without H2O Cuff Proximal Longer



Neonatal: REF 361

Pediatric: REF 371

| LENGTHS | ORDER # |
|----------|-----------|
| 2.5 x 30 | 361 - 2.5 |
| 2.5 x 38 | 371 - 2.5 |
| 3.0 x 32 | 361 - 3.0 |
| 3.0 x 39 | 371 - 3.0 |
| 3.5 x 34 | 361 - 3.5 |
| 3.5 x 40 | 371 - 3.5 |
| 4.0 x 36 | 361 - 4.0 |
| 4.0 x 41 | 371 - 4.0 |
| 4.5 x 42 | 371 - 4.5 |
| 5.0 x 44 | 371 - 5.0 |
| 5.5 x 46 | 371 - 5.5 |

Standard, with H2O Cuff



Neonatal: REF 362

Pediatric: REF 372

| LENGTHS | ORDER # |
|----------|-----------|
| 2.5 x 30 | 362 - 2.5 |
| 2.5 x 38 | 372 - 2.5 |
| 3.0 x 32 | 362 - 3.0 |
| 3.0 x 39 | 372 - 3.0 |
| 3.5 x 34 | 362 - 3.5 |
| 3.5 x 40 | 372 - 3.5 |
| 4.0 x 36 | 362 - 4.0 |
| 4.0 x 41 | 372 - 4.0 |
| 4.5 x 42 | 372 - 4.5 |
| 5.0 x 44 | 372 - 5.0 |
| 5.5 x 46 | 372 - 5.5 |

Standard, with H2O Cuff Proximal Longer



Neonatal: REF 363

Pediatric: REF 373

| LENGTHS | ORDER # |
|----------|-----------|
| 2.5 x 30 | 363 - 2.5 |
| 2.5 x 38 | 373 - 2.5 |
| 3.0 x 32 | 363 - 3.0 |
| 3.0 x 39 | 373 - 3.0 |
| 3.5 x 34 | 363 - 3.5 |
| 3.5 x 40 | 373 - 3.5 |
| 4.0 x 36 | 363 - 4.0 |
| 4.0 x 41 | 373 - 4.0 |
| 4.5 x 42 | 373 - 4.5 |
| 5.0 x 44 | 373 - 5.0 |
| 5.5 x 46 | 373 - 5.5 |

Custom Lengths (mm)

TO ORDER, SPECIFY: REF # + SIZE + LENGTH (e.g. 370 - 3.5 - C38)

Custom, without H2O Cuff



Custom: REF 370

Custom, without H2O Cuff Proximal Longer



Custom: REF 371

Custom, with H2O Cuff



Custom: REF 372

Custom, with H2O Cuff Proximal Longer



Custom: REF 373

SIZE 2.5
OD 4.4 mm

Length (mm)

| |
|--------------|
| C 26 |
| C 28 |
| 30 Std. Neo. |
| C 32 |
| C 34 |
| C 36 |
| 38 Std. Ped. |
| C 40 |
| C 42 |
| C 44 |

SIZE 3.0
OD 4.9 mm

Length (mm)

| |
|--------------|
| C 24 |
| C 26 |
| C 28 |
| C 30 |
| 32 Std. Neo. |
| C 34 |
| C 37 |
| 39 Std. Ped. |
| C 41 |
| C 43 |
| C 46 |
| C 50 |
| C 55 |

SIZE 3.5
OD 5.4 mm

Length (mm)

| |
|--------------|
| C 28 |
| C 30 |
| C 32 |
| 34 Std. Neo. |
| C 36 |
| C 38 |
| 40 Std. Ped. |
| C 42 |
| C 45 |
| C 50 |
| C 55 |
| C 60 |
| C 65 |

SIZE 4.0
OD 6.0 mm

Length (mm)

| |
|--------------|
| C 32 |
| C 34 |
| 36 Std. Neo. |
| C 38 |
| C 40 |
| 41 Std. Ped. |
| C 43 |
| C 45 |
| C 48 |
| C 51 |
| C 55 |
| C 60 |
| C 65 |

SIZE 4.5
OD 6.6 mm

Length (mm)

| |
|--------------|
| C 38 |
| C 40 |
| 42 Std. Ped. |
| C 45 |
| C 48 |
| C 51 |
| C 55 |
| C 60 |
| C 65 |

SIZE 5.0
OD 7.3 mm

Length (mm)

| |
|--------------|
| C 39 |
| C 42 |
| 44 Std. Ped. |
| C 46 |
| C 49 |
| C 53 |
| C 58 |
| C 63 |
| C 68 |

SIZE 5.5
OD 8.4 mm

Length (mm)

| |
|--------------|
| C 42 |
| C 44 |
| 46 Std. Ped. |
| C 48 |
| C 51 |
| C 55 |
| C 60 |
| C 65 |
| C 70 |